

## HRP 标记的山羊抗小鼠 IgG

中文名称： HRP 标记的山羊抗小鼠 IgG

英文名称： HRP-conjugated Goat Anti-Mouse IgG

相关类别： HRP 标记二抗

储 存： 冷冻 (-20℃)

标 记 物： HRP

### 技术规格

<b>Background</b>	Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non- immunoglobulin serum proteins. The antibody may cross-react with from other species.
<b>Applications</b>	1/ 2000-50000 for ELISA 1/ 200-2000 for Immunohisto/cytochemistry 1/ 5000-100000 for Western blotting with ECL substrates