

HRP 标记的兔抗水貂 IgG

中文名称: HRP 标记的兔抗水貂 IgG

英文名称: HRP-conjugated Rabbit anti-Mink IgG

相关类别: HRP 标记二抗

储 存: 冷冻 (-20℃)

标 记 物: HRP

技术规格

Background	Immunoglobulin G (IgG) ,is one of the most abundant proteins in serum with normal levels between 9.5-12.5 mg/ml in adult blood. Ig G is important for our defence against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by Somatic recombination and the number of specificities in an individual at a given time point is estimated to be 10 ¹¹ variants.
Applications	1/1000-5000 for ELISA 1/1000-5000 for Western blotting 1/200-500 for Immunohistochemistry 1/200-500 for Immunocytochemistry