

兔抗 NEUROG2 多克隆抗体

中文名称：兔抗 NEUROG2 多克隆抗体

英文名称：Anti-NEUROG2 rabbit polyclonal antibody

别名：NGN2, Atoh4, ngn-2, Math4A, bHLHa8

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NEUROG2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a neural-specific basic helix-loop-helix (bHLH) transcription factor that can specify a neuronal fate on ectodermal cells and is expressed in neural progenitor cells within the developing central and peripheral nervous systems. The protein product of this gene also plays a role in the differentiation and survival of midbrain dopaminergic neurons.
Applications:	ELISA, IHC
Name of antibody:	NEUROG2
Immunogen:	Synthetic peptide of human NEUROG2
Full name:	neurogenin 2

Synonyms :	NGN2, Atoh4, ngn-2, Math4A, bHLHa8
SwissProt:	Q9H2A3
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human colon cancer and Human liver cancer
IHC Recommend dilution:	15-50

