

兔抗 IFRD1 多克隆抗体

中文名称： 兔抗 IFRD1 多克隆抗体

英文名称： Anti-IFRD1 rabbit polyclonal antibody

别名： PC4; TIS7

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： IFRD1

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is an immediate early gene that encodes a protein related to interferon-gamma. This protein may function as a transcriptional co-activator/repressor that controls the growth and differentiation of specific cell types during embryonic development and tissue regeneration. Mutations in this gene are associated with sensory/motor neuropathy with ataxia. This gene may also be involved in modulating the pathogenesis of cystic fibrosis lung disease. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, IHC

Name of antibody:	IFRD1
Immunogen:	Fusion protein of human IFRD1
Full name:	interferon-related developmental regulator 1
Synonyms :	PC4; TIS7
SwissProt:	O00458
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human prostate cancer and human brain
IHC Recommend dilution:	50-200

