

## 兔抗 NEUROG3 多克隆抗体

中文名称：兔抗 NEUROG3 多克隆抗体

英文名称： Anti-NEUROG3 rabbit polyclonal antibody

别名： neurogenin 3; ngn3; Atoh5; NGN-3; Math4B; bHLHa7

抗原： NEUROG3

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

|                          |   |
|--------------------------|---|
| <b>Background:</b>       | The protein encoded by this gene is a basic helix-loop-helix (bHLH) transcription factor involved in neurogenesis. The encoded protein likely acts as a heterodimer with another bHLH protein. Defects in this gene are a cause of congenital malabsorptive diarrhea 4 (DIAR4). |
| <b>Applications:</b>     | ELISA, IHC  |
| <b>Name of antibody:</b> | NEUROG3   |
| <b>Immunogen:</b>        | Fusion protein of human NEUROG3   |
| <b>Full name:</b>        | neurogenin 3  |
| <b>Synonyms:</b>         | ngn3; Atoh5; NGN-3; Math4B; bHLHa7  |

|                                    |                        |
|------------------------------------|------------------------|
| <b>SwissProt:</b>                  | Q9Y4Z2                 |
| <b>ELISA Recommended dilution:</b> | 5000-10000             |
| <b>IHC positive control:</b>       | Human esophagus cancer |
| <b>IHC Recommend dilution:</b>     | 50-200                 |

