

## 兔抗 NF1 多克隆抗体

中文名称: 兔抗 NF1 多克隆抗体

英文名称: Anti-NF1 rabbit polyclonal antibody

级 别: WSS; NFNS; VRNF

储 存: 冷冻 (-20℃)

抗 原: NF1

宿 主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene product appears to function as a negative regulator of the ras signal transduction pathway. Mutations in this gene have been linked to neurofibromatosis type 1, juvenile myelomonocytic leukemia and Watson syndrome. The mRNA for this gene is subject to RNA editing (CGA>UGA->Arg1306Term) resulting in premature translation termination. Alternatively spliced transcript variants encoding different isoforms have also been described for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NF1

<b>Immunogen:</b>	Synthetic peptide of human NF1
<b>Full name:</b>	neurofibromin 1
<b>Synonyms:</b>	WSS; NFNS; VRNF
<b>SwissProt:</b>	P21359
<b>ELISA Recommended dilution:</b>	500-1000
<b>IHC positive control:</b>	Human liver cancer and human gastric cancer
<b>IHC Recommend dilution:</b>	40-250



