

HRP 标记山羊抗兔/小鼠 IgG（免疫组化，即用型）

中文名称：HRP 标记山羊抗兔/小鼠 IgG（免疫组化，即用型）

英文名称：HRP-conjugated Goat Anti-Rabbit/mouse IgG for IHC(ready to use)

相关类别：HRP 标记二抗

储 存：冷藏：（2-8℃）

标 记 物：HRP

适用样本：

Specimens: Formalin-fixed, paraffin-embedded tissue sections as well as frozen tissue sections and cell smears.

技术规格

Applications:	This product is well suited for a two-step immunohistochemistry procedure. And it is provided ready to use. This product reacts with rabbit and mouse immunoglobulins of all classes. Further dilution of this reagent or alteration of the incubation time or temperature may give erroneous results.
---------------	--