

兔抗 UBA2 多克隆抗体

中文名称：兔抗 UBA2 多克隆抗体

英文名称： Anti-UBA2 rabbit polyclonal antibody

别名： ARX; SAE2; HRIHFB2115

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： UBA2

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Posttranslational modification of proteins by the addition of the small protein SUMO (see SUMO1; MIM 601912), or sumoylation, regulates protein structure and intracellular localization. SAE1 (MIM 613294) and UBA2 form a heterodimer that functions as a SUMO-activating enzyme for the sumoylation of proteins (Okuma et al., 1999 [PubMed 9920803]).
Applications:	WB, IHC
Name of antibody:	UBA2
Immunogen:	Synthesized peptide derived from C-terminal o

	f human Uba2.
Full name:	ubiquitin-like modifier activating enzyme 2
Synonyms :	ARX; SAE2; HRIHFB2115
SwissProt:	Q9UBT2
IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	71 kDa
WB Positive control:	293 cells lysates treated with UV (5mins)
WB Recommended dilution:	500-3000



