

HRP 标记兔抗大鼠 IgG

中文名称: HRP 标记兔抗大鼠 IgG

英文名称: HRP-conjugated Rabbit anti-rat IgG

相关类别: HRP 标记二抗 储存 冷冻 (-20℃)

标记物: HRP

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule rat IgG. It also reacts with the light chains of other rat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
Applications	1/2000-50000 for ELISA 1/200-2000 for Immunohisto/cytochemistry 1/5000-100000 for Western blotting with ECL substrates