

## HRP 标记兔抗狗 IgG

中文名称: HRP 标记兔抗狗 IgG

英文名称: HRP-conjugated Rabbit Anti-Dog IgG

相关类别: HRP 标记二抗

储 存: 冷冻(-20℃)

标 记 物: HRP

## 技术规格

Background:	Immunoglobulin G (IgG) is one of the most abundant proteins in seru m with normal levels between 8-17 mg/mL in adult blood. IgG is imp ortant for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune respon se. The IgG molecule has two separate functions, to bind to the path ogen that elicited the response and to recruit other cells and molecul es to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.
Applications:	1/2000-10000 for ELISA 1/1000-5000 for Western blotting 1/200-1000 for Immunohistochemistry/cytochemistry