

## 兔抗 UBA5 多克隆抗体

中文名称: 兔抗 UBA5 多克隆抗体

英文名称: Anti-UBA5 rabbit polyclonal antibody

别名: THIFP1; UBE1DC1

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: UBA5

反应种属: Human, Mouse, Rat

标记物: Unconjugate

相关类别: 一抗

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the E1-like ubiquitin-activating enzyme family. This protein activates ubiquitin-fold modifier 1, a ubiquitin-like post-translational modifier protein, via the formation of a high-energy thioester bond. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. A pseudogene located on chromosome 1 has also been identified.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UBA5
<b>Immunogen:</b>	Fusion protein of human UBA5

<b>Full name:</b>	ubiquitin-like modifier activating enzyme 5
<b>Synonyms:</b>	THIFP1; UBE1DC1
<b>SwissProt:</b>	Q9GZZ9
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and human colorectal cancer
<b>IHC Recommend dilution:</b>	25-100



