

HRP 标记兔抗人 IgG

中文名称： HRP 标记兔抗人 IgG

英文名称： HRP-conjugated Rabbit anti-human IgG

相关类别： HRP 标记二抗

储 存： 冷冻 (-20℃)

标 记 物： HRP

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule human IgG. It also reacts with the light chains of other human immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
Applications	1/1000-20000 for ELISA 1/100-1000 for Immunohisto/cytochemistry 1/2000-40000 for Western blotting with ECL substrates