

HRP 标记兔抗山羊 IgG

中文名称： HRP 标记兔抗山羊 IgG

英文名称： HRP-conjugated Rabbit anti-Goat IgG

相关类别： HRP 标记二抗

储 存： 冷冻（-20℃）

标 记 物： HRP

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule goat IgG. It also reacts with the light chains of other goat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.
Applications	1/2000-50000 for ELISA 1/200-2000 for Immunohisto/cytochemistry 1/5000-100000 for Western blotting with ECL substrates