

HRP 标记小鼠抗兔 IgG

中文名称: HRP 标记小鼠抗兔 IgG

英文名称: HRP-conjugated Mouse anti-rabbit IgG

相关类别: HRP 标记二抗

储 存: 冷冻(-20℃)

标 记 物: HRP

技术规格

Background	Antibody Specificity: Based on immunoelectrophoresis and/or ELIS A, the antibody reacts with whole molecule rabbit IgG. It also reacts with the light chains of other rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody exhibits inherent minimal cross-reaction to mouse serum proteins and has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with human, goat, mouse, and sheep serum proteins. The antibody may cross-react with immunoglobulins from other species.
Applications	1/2000-50000 for ELISA 1/200-2000 for Immunohisto/cytochemistr y 1/5000-100000 for Western blotting with ECL substrates