

## 兔抗 UBA7 多克隆抗体

中文名称: 兔抗 UBA7 多克隆抗体

英文名称: Anti-UBA7 rabbit polyclonal antibody

别名: D8; UBE2; UBA1B; UBE1L

抗原: UBA7

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

**Background:**

The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E1 ubiquitin-activating enzyme family. The encoded enzyme is a retinoid target that triggers promyelocytic leukemia (PML)/retinoic acid receptor alpha (RARalpha) degradation and apop

	tosis in acute promyelocytic leukemia, where it is involved in the conjugation of the ubiquitin-like interferon-stimulated gene 15 protein.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	UBA7
<b>Immunogen:</b>	Fusion protein of human UBA7
<b>Full name:</b>	ubiquitin-like modifier activating enzyme 7
<b>Synonyms:</b>	D8; UBE2; UBA1B; UBE1L
<b>SwissProt:</b>	P41226
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and human esophagus cancer
<b>IHC Recommend dilution:</b>	30-150



