

PE 标记的小鼠抗山羊 IgG

- 中文名称: PE 标记的小鼠抗山羊 lgG
- 英文名称: PE-conjugated Mouse anti-goat IgG
- 相关类别: PE 标记二抗
- 储 存: 冷冻 (-20℃)
- 标记物: PE

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins in seru m with normal levels between 8-17 mg/ml in adult blood. IgG is impo rtant for our defense against microorganisms and the molecules are pr oduced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by so matic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.
Applications	1/20-100 for FCM 1/20-200 for Immunohisto/cytochemistry