

PE 标记山羊抗兔 IgG

中文名称: PE 标记山羊抗兔 IgG

英文名称: PE-conjugated Goat anti-Rabbit IgG

相关类别: PE 标记二抗

储 存: 冷冻(-20℃)

标 记 物: PE

技术规格

| Background | Immunoglobulin G (IgG) is one of the most abundant proteins in |
|--------------|--|
| | serum with normal levels between 8-17 mg/mL in adult blood. Ig |
| | G is important for our defense against microorganisms and the m |
| | olecules are produced by B lymphocytes as a part of our adaptive |
| | immune response. The IgG molecule has two separate functions, t |
| | o bind to the pathogen that elicited the response and to recruit |
| | other cells and molecules to destroy the antigen. The variability o |
| | f the IgG pool is generated by somatic recombination and the nu |
| | mber of specificities in an individual at a given time point is esti |
| | mated to be 1011 variants. |
| Applications | 1/50-200 for most application |

