

兔抗 IGF1R 多克隆抗体

- 中文名称： 兔抗 IGF1R 多克隆抗体
- 英文名称： Anti-IGF1R (Ab-1165/1166) rabbit polyclonal antibody
- 别名： IGF1R; CD221; IGFIR; JTK13
- 相关类别： 一抗
- 储存： 冷冻（-20℃） 避光
- 宿主： Rabbit
- 抗原： IGF1R
- 反应种属： Human, Mouse, Rat
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

技术规格

Background:	This receptor binds insulin-like growth factor 1 (IGF1) with a high affinity and IGF2 with a lower affinity. It has a tyrosine-protein kinase activity, which is necessary for the activation of the IGF1-stimulated downstream signaling cascade. When present in a hybrid receptor with INSR, binds IGF1.
Applications:	WB
Name of antibody:	IGF1R
Immunogen:	Synthetic peptide of human IGF1R (Ab-1165/1166)
Full name:	insulin-like growth factor 1 receptor
Synonyms :	IGFR; CD221; IGFIR; JTK13

SwissProt:	P08069
WB Predicted band size:	155 kDa
WB Positive control:	293 cells
WB Recommended dilution:	500-1000

