

IFNW1 抗原（重组蛋白）

中文名称： IFNW1 抗原（重组蛋白）

英文名称： IFNW1 Antigen (Recombinant Protein)

别名： interferon omega 1

相关类别： 抗原

储存： 冷冻（-20℃）

概述

Fusion protein corresponding to a region derived from 1-195 amino acids of human IFNW1

技术规格

Full name:	interferon omega 1
Swissprot:	P05000
Gene Accession:	BC117290
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is an interferon and possesses a antiviral activity. The encoded protein binds to the interferon alpha/beta receptor but not to the interferon gamma receptor. This intro nless gene has several pseudogenes spread throughout the genom e.

mlbio 梅联生物
Good elisakit producers