

## 兔抗 NFE2L1 多克隆抗体

中文名称：兔抗 NFE2L1 多克隆抗体

英文名称： Anti-NFE2L1 rabbit polyclonal antibody

别名： NRF1; TCF11; LCR-F1

储存： 冷冻（-20℃）

抗原： NFE2L1

宿主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that is involved in globin gene expression in erythrocytes. Confusion has occurred in bibliographic databases due to the shared symbol of NRF1 for this gene, NFE2L1, and for "nuclear respiratory factor 1" which has an official symbol of NRF1.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	NFE2L1
<b>Immunogen:</b>	Fusion protein of human NFE2L1
<b>Full name:</b>	nuclear factor, erythroid 2-like 1
<b>Synonyms :</b>	NRF1; TCF11; LCR-F1
<b>SwissProt:</b>	Q14494

<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human thyroid cancer and human brain
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	85 kDa
<b>WB Positive control:</b>	Raji, hepG2 and hela cells
<b>WB Recommended dilution:</b>	200-1000



