

金标记山羊抗兔 IgG

中文名称: 金标记山羊抗兔 IgG

英文名称: Gold-conjugated Goat anti-rabbit IgG

相关类别: 其他二抗

储 存: 冷冻(-20℃)

标记物: Gold

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant p roteins in serum with normal levels between 8-17 mg/m I in adult blood. IgG is important for our defense again st microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the response and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by somatic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.
Applications	1/20-100 for immunoelectron microscopy(IEM), immunog old staining(IGS), gold immunochromatography assay(GI CA)
Conjugate	Gold (15 nm)

