

兔抗 REL 多克隆抗体

中文名称：兔抗 REL 多克隆抗体

英文名称： Anti-REL rabbit polyclonal antibody

别名： C-Rel

相关类别： 一抗

抗原： REL

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The REL gene encodes c-Rel, a transcription factor that is a member of the Rel/NFKB family, which also includes RELA (MIM 164014), RELB (604758), NFKB1 (MIM 164011), and NFKB2 (MIM 164012). These proteins are related through a highly conserved N-terminal region termed the 'Rel domain,' which is responsible for DNA binding, dimerization, nuclear localization, and binding to the NFKB inhibitor
Applications:	ELISA, WB, IHC
Name of antibody:	REL
Immunogen:	Fusion protein of human REL

Full name:	v-rel avian reticuloendotheliosis viral oncogene homolog
Synonyms :	C-Rel
SwissProt:	Q04864
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer and Human brain
IHC Recommend dilution:	100-300
WB Predicted band size:	69 kDa
WB Positive control:	Hela and 293T cells
WB Recommended dilution:	1000-5000



