

兔抗 REM1 多克隆抗体

中文名称：兔抗 REM1 多克隆抗体

英文名称：Anti-REM1 rabbit polyclonal antibody

别名：RRAD and GEM like GTPase 1 GES; GD:REM

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：REM1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a GTPase and member of the RAS-like GTP-binding protein family. The encoded protein is expressed in endothelial cells, where it promotes reorganization of the actin cytoskeleton and morphological changes in the cells.
Applications:	ELISA, IHC
Name of antibody:	REM1
Immunogen:	Synthetic peptide of human REM1
Full name:	RRAD and GEM like GTPase 1
Synonyms:	GES; GD:REM

SwissProt:	O75628
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer;
IHC Recommend dilution:	20-100



