

IL22RA2 抗原（重组蛋白）

中文名称：IL22RA2 抗原（重组蛋白）

英文名称： IL22RA2 Antigen (Recombinant Protein)

别名： interleukin 22 receptor subunit alpha 2; CRF2X; CRF2-10; CRF2-S1; IL-22BP; IL-22RA2; ZCYTOR16; IL-22R-alpha-2

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 22-231 amino acids of human IL22RA2

技术规格

Full name:	interleukin 22 receptor subunit alpha 2
Synonyms:	CRF2X; CRF2-10; CRF2-S1; IL-22BP; IL-22RA2; ZCYTOR16; IL-22R-alpha-2
Swissprot:	Q969J5
Gene Accession:	BC125167
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the class II cytokine receptor family. The encoded soluble protein specifically binds to and inhibits interleukin 22 activity by blocking the interaction of interleukin 22 with its cell surface receptor. The encoded protein may be important in the regulation of inflammatory response, and has been implicated

d in the regulation of tumorigenesis in the colon. Alternative splicing results in multiple transcript variants encoding distinct isoforms.