

兔抗 REM2 多克隆抗体

中文名称：兔抗 REM2 多克隆抗体

英文名称：Anti-REM2 rabbit polyclonal antibody

别名：AW411893

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：REM2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	REM2, a Ras-like GTPase, may be involved in increasing endothelial cell proliferation and the induction of angiogenesis. REM2 belongs to the small GTPase superfamily, RGK family. It binds GTP saturably and exhibits a low intrinsic rate of GTP hydrolysis.
Applications:	ELISA, WB, IHC, IF
Name of antibody:	REM2
Immunogen:	Synthetic peptide of human REM2
Full name:	rad and gem related GTP binding protein 2
Synonyms :	AW411893
SwissProt:	Q8IYK8

ELISA Recommended dilution:	1000-2000
IHC positive control:	Human spleen
IHC Recommend dilution:	200-400
WB Predicted band size:	37 kDa
WB Positive control:	NCCIT cells
WB Recommended dilution:	200-1000
IF positive control:	F9 cells
IF Recommended dilution:	20-50



