

链霉亲和素标记的兔抗人 IgG

- 中文名称: 链霉亲和素标记的兔抗人 lgG
- 英文名称: Streptavidin-conjugated Rabbit anti-human IgG
- 相关类别: 其他二抗
- 储 存: 冷冻 (-20℃)
- 标 记 物: Streptavidin

技术规格

Background	Immunoglobulin G (IgG), is one of the most abundant proteins in serum with normal levels between 9.5-12.5 mg/ml in adult blood. IgG is impor tant for our defense against microorganisms and the molecules are pro duced by B lymphocytes as a part of our adaptive immune response. T he IgG molecule has two separate functions; to bind to the pathogen t hat elicited the response and to recruit other cells and molecules to de stroy the antigen. The variability of the IgG pool is generated by somati c recombination and the number of specificities in an individual at a giv en time point is estimated to be 1011 variants.
Applications	1/500-5000 for ELISA 1/50-50 for Immunohistochemistry/cytochemistry