

## 兔抗 REPS1 多克隆抗体

中文名称: 兔抗 REPS1 多克隆抗体

英文名称: Anti-REPS1 rabbit polyclonal antibody

别名: RALBP1

抗原: REPS1

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human Rat

标记物: Unconjugate

相关类别: 一抗

克隆类型: rabbit polyclonal

### 技术规格

#### Background:

This gene encodes a signaling adaptor protein with two E H domains that interacts with proteins that participate in signaling, endocytosis and cytoskeletal changes. The encoded protein has been found in association with intersectin 1 and Src homology 3-domain growth factor receptor-bound 2-like (endophilin) interacting protein 1 when intersectin 1 was isolated from clathrin-coated pits. The encoded protein has also been shown to interact with amphiphysin, a cytoplasmic protein at the surface of synaptic vesicles. Alternative splicing results in multiple transcript variants encoding different isoforms.

<b>Applications:</b>	WB
<b>Name of antibody:</b>	REPS1
<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human REPS1.
<b>Full name:</b>	RALBP1 associated Eps domain containing 1
<b>Synonyms :</b>	RALBP1
<b>SwissProt:</b>	Q96D71
<b>WB Predicted band size:</b>	87 kDa
<b>WB Positive control:</b>	Rat heart tissue lysate
<b>WB Recommended dilution:</b>	500-3000

