

## IL31RA 抗原（重组蛋白）

中文名称： IL31RA 抗原（重组蛋白）

英文名称： IL31RA Antigen (Recombinant Protein)

别名： interleukin 31 receptor A; CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to C terminal 190 amino acids of human IL31RA

### 技术规格

<b>Full name:</b>	interleukin 31 receptor A
<b>Synonyms:</b>	CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384
<b>Swissprot:</b>	Q8NI17
<b>Gene Accession:</b>	BC110490
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene belongs to the type I cytokine receptor family. This receptor, with homology to gp130, is expressed on monocytes, and is involved in IL-31 signaling via activation of STAT-3 and STAT-5. It functions either as a monomer, or as part of a receptor complex with oncostatin M receptor (OSMR). Several alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

**mlbio** 梅联生物  
Good elisakit producers