

## 驴抗山羊 IgG

中文名称： 驴抗山羊 IgG

英文名称： Donkey anti-Goat IgG

相关类别： 其他二抗

储 存： 冷冻 (-20℃)

标 记 物： Unconjugate

### 技术规格

#### Background

Antibody Specificity: Based on immunoelectrophoresis, the antibody reacts with whole molecule goat IgG. It also reacts with the light chains of other goat immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with immunoglobulins from other species.