

山羊抗牛 IgG (亲和纯化)

中文名称: 山羊抗牛 IgG (亲和纯化)

英文名称: Goat anti-Bovine IgG (AffiniPure)

相关类别: 其他二抗

储 存: 冷冻(-20℃)

技术规格

Background	Immunoglobulin G (IgG) is one of the most abundant proteins in serum with normal levels between 9.5-12.5 mg/ml in adult blood. IgG is important for our defense against microorganism s and the molecules are produced by B lymphocytes as a part of our adaptive immune response. The IgG molecule has two separate functions; to bind to the pathogen that elicited the r esponse and to recruit other cells and molecules to destroy the antigen. The variability of the IgG pool is generated by som atic recombination and the number of specificities in an individual at a given time point is estimated to be 1011 variants.
Applications	1/1000-5000 for ELISA 1/1000-5000 for Western blotting 1/200 -500 for Immunohistochemistry/Immunocytochemistry