

## 兔抗 NFKBIA(Ab-32/36) 多克隆抗体

中文名称：兔抗 NFKBIA(Ab-32/36) 多克隆抗体

英文名称： Anti-NFKBIA(Ab-32/36) rabbit polyclonal antibody

别名： IKBA; MAD-3; NFKBI

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： NFKBIA(Ab-32/36)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the NF-kappa-B inhibitor family, which contain multiple ankrin repeat domains. The encoded protein interacts with REL dimers to inhibit NF-kappa-B/REL complexes which are involved in inflammatory responses. The encoded protein moves between the cytoplasm and the nucleus via a nuclear localization signal and CRM1-mediated nuclear export. Mutations in this gene have been found in ectodermal dysplasia anhidrotic with T-cell immunodeficiency autosomal dominant disease.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	NFKBIA(Ab-32/36)

<b>Immunogen:</b>	Synthesized non-phosphopeptide derived from human I $\kappa$ B- $\alpha$ around the phosphorylation site of Serine 32/36.
<b>Full name:</b>	nuclear factor of kappa light polypeptide gene enhancer in B-cells inhibitor, alpha
<b>Synonyms :</b>	IKBA; MAD-3; NFKBI
<b>SwissProt:</b>	P25963
<b>IHC positive control:</b>	Human breast carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	36 kDa
<b>WB Positive control:</b>	MCF7 cells treated with TNF- $\alpha$ lysate
<b>WB Recommended dilution:</b>	500-3000



