

兔抗 RETN 多克隆抗体

中文名称：兔抗 RETN 多克隆抗体

英文名称： Anti-RETN rabbit polyclonal antibody

别名： ADSF, RSTN, XCP1, FIZZ3, RETN1

抗原： RETN

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene belongs to the family defined by the mouse resistin-like genes. The characteristic feature of this family is the C-terminal stretch of 10 cys residues with identical spacing. The mouse homolog of this protein is secreted by adipocytes, and may be the hormone potentially linking obesity to type II diabetes. Alternatively spliced transcript variants encoding the same protein have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	RETN
Immunogen:	Fusion protein of human RETN

Full name:	resistin
Synonyms :	ADSF, RSTN, XCP1, FIZZ3, RETN1
SwissProt:	Q9HD89
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human ovarian cancer and Human lung cancer
IHC Recommend dilution:	25-100

