

兔抗 NFKBIZ 多克隆抗体

中文名称：兔抗 NFKBIZ 多克隆抗体

英文名称： Anti-NFKBIZ rabbit polyclonal antibody

别名： IKBZ; INAP; MAIL

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NFKBIZ

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene is a member of the ankyrin-repeat family and is induced by lipopolysaccharide (LPS). The C-terminal portion of the encoded product which contains the ankyrin repeats, shares high sequence similarity with the I kappa B family of proteins. The latter are known to play a role in inflammatory responses to LPS by their interaction with NF-B proteins through ankyrin-repeat domains. Studies in mouse indicate that this gene product is one of the nuclear I kappa B proteins and an activator of IL-6 production. Two transcript variants encoding different isoforms have been found for this gene.

Applications:	ELISA, WB
Name of antibody:	NFKBIZ
Immunogen:	Synthetic peptide of human NFKBIZ
Full name:	NFKB inhibitor zeta
Synonyms:	IKBZ; INAP; MAIL
SwissProt:	Q9BYH8
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	78 kDa
WB Positive control:	293T cell lysate
WB Recommended dilution:	500-2000

