

兔抗 NFYA 多克隆抗体

中文名称：兔抗 NFYA 多克隆抗体

英文名称： Anti-NFYA rabbit polyclonal antibody

别名： HAP2; CBF-A; CBF-B; NF-YA

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： NFYA

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is one subunit of a trimeric complex, forming a highly conserved transcription factor that binds to CCAAT motifs in the promoter regions in a variety of genes. Subunit A associates with a tight dimer composed of the B and C subunits, resulting in a trimer that binds to DNA with high specificity and affinity. The sequence specific interactions of the complex are made by the A subunit, suggesting a role as the regulatory subunit.
Applications:	WB
Name of antibody:	NFYA
Immunogen:	Synthesized peptide derived from internal of human NFYA.

Full name:	nuclear transcription factor Y, alpha
Synonyms :	HAP2; CBF-A; CBF-B; NF-YA
SwissProt:	P23511
WB Predicted band size:	37 kDa
WB Positive control:	293 cells, HepG2 cells and COLO205 cells lysates
WB Recommended dilution:	500-3000

