

兔抗 NGB 多克隆抗体

中文名称：兔抗 NGB 多克隆抗体

英文名称： Anti-NGB rabbit polyclonal antibody

抗 原： NGB

储 存： 冷冻（-20℃）

宿 主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an oxygen-binding protein that is distantly related to members of the globin gene family. It is highly conserved among other vertebrates. It is expressed in the central and peripheral nervous system where it may be involved in increasing oxygen availability and providing protection under hypoxic/ischemic conditions.
Applications:	ELISA, IHC
Name of antibody:	NGB
Immunogen:	Fusion protein of human NGB
Full name:	neuroglobin
SwissProt:	Q9NPG2
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer

IHC Recommend dilution: 25-100

