

## 兔抗 NHLRC2 多克隆抗体

中文名称：兔抗 NHLRC2 多克隆抗体

英文名称：Anti-NHLRC2 rabbit polyclonal antibody

别名：NHL repeat containing 2; FINCA

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NHLRC2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Required for normal embryonic development.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NHLRC2
<b>Immunogen:</b>	Fusion protein of human NHLRC2
<b>Full name:</b>	NHL repeat containing 2
<b>Synonyms:</b>	FINCA
<b>SwissProt:</b>	Q8NBF2
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	150-300

