

## 兔抗 NIF3L1 多克隆抗体

中文名称：兔抗 NIF3L1 多克隆抗体

英文名称：Anti-NIF3L1 rabbit polyclonal antibody

别名：CALS-7; MDS015; ALS2CR1

相关类别：一抗

抗原：NIF3L1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

<b>Background:</b>	May function as a transcriptional corepressor through its interaction with COPS2, negatively regulating the expression of genes involved in neuronal differentiation.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NIF3L1
<b>Immunogen:</b>	Fusion protein of human NIF3L1
<b>Full name:</b>	NGG1 interacting factor 3 like 1
<b>Synonyms:</b>	CALS-7; MDS015; ALS2CR1
<b>SwissProt:</b>	Q9GZT8
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human gastric cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	50-200

