

## 兔抗 NIPA1 多克隆抗体

中文名称：兔抗 NIPA1 多克隆抗体

英文名称：Anti-NIPA1 rabbit polyclonal antibody

别名：FSP3; SPG6

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NIPA1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a magnesium transporter that associates with early endosomes and the cell surface in a variety of neuronal and epithelial cells. This protein may play a role in nervous system development and maintenance. Multiple transcript variants encoding different isoforms have been found for this gene. Mutations in this gene have been associated with autosomal dominant spastic paraplegia 6.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NIPA1
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	non imprinted in Prader-Willi/Angelman syndrome 1
<b>Synonyms :</b>	FSP3; SPG6

<b>SwissProt:</b>	Q7RTP0
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	25-100

kDa

130 —

95 —

72 —

55 —

36 —

28 —

17 —



