

## 兔抗 NIPSNAP1 多克隆抗体

中文名称：兔抗 NIPSNAP1 多克隆抗体

英文名称：Anti-NIPSNAP1 rabbit polyclonal antibody

别名：nipsnap homolog 1

相关类别：一抗

抗原：NIPSNAP2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

<b>Background:</b>	This gene encodes a member of the NipSnap family of proteins that may be involved in vesicular transport. The encoded protein is localized to mitochondria and plays a role in oxidative phosphorylation. A pseudogene of this gene is located on the long arm of chromosome 2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	NIPSNAP2
<b>Immunogen:</b>	Synthetic peptide of human NIPSNAP2
<b>Full name:</b>	nipsnap homolog 2
<b>Synonyms:</b>	GBAS

<b>SwissProt:</b>	O75323
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	34 kDa
<b>WB Positive control:</b>	Mouse liver tissue, Human fetal liver tissue, HepG2 cell lysates
<b>WB Recommended dilution:</b>	200-1000



