

## 兔抗 NIPSNAP3A 多克隆抗体

中文名称: 兔抗 NIPSNAP3A 多克隆抗体

英文名称: Anti-NIPSNAP3A rabbit polyclonal antibody

别名: nipsnap homolog 3A; TASSC; HSPC299; NIPSNAP4

储存: 冷冻 (-20℃)

抗原: NIPSNAP3A

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	NIPSNAP3A belongs to a family of proteins with putative roles in vesicular transport (Buechler et al., 2004 [PubMed 15177564]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NIPSNAP3A
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	nipsnap homolog 3A
<b>Synonyms:</b>	TASSC; HSPC299; NIPSNAP4
<b>SwissProt:</b>	Q9UFN0
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	20-100

