

兔抗 NKX2-1 多克隆抗体

中文名称：兔抗 NKX2-1 多克隆抗体

英文名称： Anti-NKX2-1 rabbit polyclonal antibody

别名： BCH, BHC, NK-2, TEBP, TTF1, NKX2A, TITF1, TTF-1, NKX2.1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NKX2-1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a protein initially identified as a thyroid-specific transcription factor. The encoded protein binds to the thyroglobulin promoter and regulates the expression of thyroid-specific genes but has also been shown to regulate the expression of genes involved in morphogenesis. Mutations and deletions in this gene are associated with benign hereditary chorea, choreoathetosis, congenital hypothyroidism, and neonatal respiratory distress, and may be associated with thyroid cancer. Multiple transcript variants encoding different isoforms have been found for this gene. This gene shares the symbol/alias 'TTF1' with another gene, transcription termination factor 1, which plays a role in ribosomal gene transcription.

Applications:	ELISA, WB
Name of antibody:	NKX2-1
Immunogen:	Synthetic peptide of human NKX2-1
Full name:	NK2 homeobox 1
Synonyms :	BCH, BHC, NK-2, TEBP, TTF1, NKX2A, TITF1, TTF-1, NKX2.1
SwissProt:	P43699
ELISA Recommended dilution:	1000-5000
WB Predicted band size:	42 kDa
WB Positive control:	Mouse lung tissue
WB Recommended dilution:	200-1000

