



## 兔抗 NKX2-8 多克隆抗体

中文名称： 兔抗 NKX2-8 多克隆抗体

英文名称： Anti-NKX2-8 rabbit polyclonal antibody

别 名： NKX2H; NKX2.8; Nkx2-9

相关类别： 一抗

储 存： 冷冻 (-20°C) 避光

宿 主： Rabbit

抗 原： NKX2-8

反应种属： Human

标 记 物： Unconjugate

克隆类型： Unconjugate

技术规格

<b>Background:</b>	Members of the NK-2 family of homeodomain proteins, which include Nkx-2.2, Nkx-2.5, Nkx-2.6 and Nkx-2.8, are key regulators of growth and development in several tissues, including brain, heart and pancreas. Nkx-2.8, also designated NK-2 homolog H, NKX2H or Nkx-2.9, is a nuclear protein that contains one homeobox DNA-binding domain, indicating a possible role in development.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	NKX2-8
<b>Immunogen:</b>	Synthesized peptide derived from N-terminal of human NKX2-8.
<b>Full name:</b>	NK2 homeobox 8
<b>Synonyms :</b>	NKX2H; NKX2.8; Nkx2-9

<b>SwissProt:</b>	Q8IUT7
<b>WB Predicted band size:</b>	26 kDa
<b>WB Positive control:</b>	K562 cells lysate
<b>WB Recommended dilution:</b>	500-3000

