

兔抗 NLE1 多克隆抗体

中文名称: 兔抗 NLE1 多克隆抗体

英文名称: Anti-NLE1 rabbit polyclonal antibody

别名: Nle

储存: 冷冻 (-20℃) 避光

抗原: NLE1

宿主: Rabbit

反应种属: Human Mouse Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Plays a role in regulating Notch activity. Plays a role in regulating the expression of CDKN1A and several members of the Wnt pathway, probably via its effects on Notch activity. Required during embryogenesis for inner mass cell survival (By similarity).
Applications:	WB
Name of antibody:	NLE1
Immunogen:	Synthesized peptide derived from internal of human NLE1.
Full name:	notchless homolog 1
Synonyms:	NLE
SwissProt:	Q9NVX2
WB Predicted band size:	53 kDa
WB Positive control:	COS cells lysate

WB Recommended dilution: 500-3000

