

## 兔抗 NME1 多克隆抗体

中文名称：兔抗 NME1 多克隆抗体

英文名称： Anti-NME1 rabbit polyclonal antibody

别名： NB, AWD, NBS, GAAD, NDKA, NM23, NDPKA, NDPK-A, NM23-H1

抗原： NME1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

<b>Background:</b>	<p>This gene (NME1) was identified because of its reduced mRNA transcript levels in highly metastatic cells. Nucleoside diphosphate kinase (NDK) exists as a hexamer composed of 'A' (encoded by this gene) and 'B' (encoded by NME2) isoforms. Mutations in this gene have been identified in aggressive neuroblastomas. Two transcript variants encoding different isoforms have been found for this gene. Co-transcription of this gene and the neighboring downstream gene (NME2) generates naturally-occurring transcripts (NME1-NME2), which encodes a fusion protein comprised of sequence sharing identity with each individual gene product.</p>
--------------------	---

<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	NME1
<b>Immunogen:</b>	Fusion protein of human NME1
<b>Full name:</b>	NME/NM23 nucleoside diphosphate kinase 1
<b>Synonyms :</b>	NB, AWD, NBS, GAAD, NDKA, NM23, NDPKA, NDPK-A, NM23-H1
<b>SwissProt:</b>	P15531
<b>ELISA Recommended dilution:</b>	2000-10000
<b>WB Predicted band size:</b>	17 kDa
<b>WB Positive control:</b>	A549 cells
<b>WB Recommended dilution:</b>	1000-5000

