

兔抗 NME2 多克隆抗体

中文名称：兔抗 NME2 多克隆抗体

英文名称：Anti-NME2 rabbit polyclonal antibody

别名：PUF; NDKB; NDPKB; NM23B; NDPK-B; NM23-H2

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：NME2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Nucleoside diphosphate kinase (NDK) exists as a hexamer composed of 'A' (encoded by NME1) and 'B' (encoded by this gene) isoforms. Multiple alternatively spliced transcript variants have been found for this gene. Read-through transcription from the neighboring upstream gene (NME1) generates naturally-occurring transcripts (NME1-NME2) that encode a fusion protein comprised of sequence sharing identity with each individual gene product.
Applications:	ELISA, WB, IHC
Name of antibody:	NME2

Immunogen:	Synthetic peptide of humanNME2
Full name:	NME/NM23 nucleoside diphosphate kinase 2
Synonyms :	PUF; NDKB; NDPKB; NM23B; NDPK-B; NM23-H2
SwissProt:	P22392
ELISA Recommended dilution :	2000-10000
IHC positive control:	Human thyroid cancer and Human cervical cancer
IHC Recommend dilution:	30-150
WB Predicted band size:	17 kDa
WB Positive control:	Jurkat and A549 cells, Mouse brain tissue, 293T, LNCap, NI H/3T3 and Hela cells
WB Recommended dilution:	1000-5000



