

## 兔抗 NME4 多克隆抗体

中文名称：兔抗 NME4 多克隆抗体

英文名称：Anti-NME4 rabbit polyclonal antibody

别名：NME/NM23 nucleoside diphosphate kinase 4 NDPK-D; NM23H4; nm23-H4

抗原：NME4

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

<b>Background:</b>	The nucleoside diphosphate (NDP) kinases (EC 2.7.4.6) are ubiquitous enzymes that catalyze transfer of gamma-phosphates, via a phosphohistidine intermediate, between nucleoside and dioxynucleoside tri- and diphosphates. The enzymes are products of the nm23 gene family, which includes NME4 (Milon et al., 1997 [PubMed 9099850]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NME4
<b>Immunogen:</b>	Fusion protein of human NME4
<b>Full name:</b>	NME/NM23 nucleoside diphosphate kinase 4
<b>Synonyms:</b>	NDPK-D; NM23H4; nm23-H4

<b>SwissProt:</b>	O00746
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer; Human gastric cancer



