

兔抗 NME6 多克隆抗体

中文名称：兔抗 NME6 多克隆抗体

英文名称： Anti-NME6 rabbit polyclonal antibody

别名： NDK 6; NM23-H6; IPIA-ALPHA

抗原： NME6

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Nucleoside diphosphate (NDP) kinases (EC 2.7.4.6), such as NME6, are ubiquitous enzymes that catalyze transfer of gamma-phosphates, via a phosphohistidine intermediate, between nucleoside and dioxynucleoside tri- and diphosphates. Major role in the synthesis of nucleoside triphosphates other than ATP. The ATP gamma phosphate is transferred to the NDP beta phosphate via a ping-pong mechanism, using a phosphorylated active-site intermediate. Inhibitor of p53-induced apoptosis.
Applications:	ELISA, WB
Name of antibody:	NME6

Immunogen:	Fusion protein of human NME6
Full name:	NME/NM23 nucleoside diphosphate kinase 6
Synonyms :	NDK 6; NM23-H6; IPIA-ALPHA
SwissProt:	O75414
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	21 kDa
WB Positive control:	Mouse spleen and intestinum crassum tissue
WB Recommended dilution:	200-1000

