

## 山羊抗兔 IgG

中文名称: 山羊抗兔 lgG

英文名称: Goat anti-Rabbit IgG

相关类别: 其他二抗

储 存: 冷冻(-20℃)

## 技术规格

Immunoglobulin G (IgG) is one of the most abundant proteins in seru m with normal levels between 9.5-12.5 mg/ml in adult blood. IgG is im portant for our defense against microorganisms and the molecules are produced by B lymphocytes as a part of our adaptive immune respons e. The IgG molecule has two separate functions; to bind to the pathog en that elicited the response and to recruit other cells and molecules t o destroy the antigen. The variability of the IgG pool is generated by s omatic recombination and the number of specificities in an individual a
t a given time point is estimated to be 1011 variants.  Immune colloidal gold immunochromatography