

兔抗 NMS 多克隆抗体

中文名称: 兔抗 NMS 多克隆抗体

英文名称: Anti-NMS rabbit polyclonal antibody

储 存: 冷冻 (-20℃) 避光

抗 原: NMS

宿 主: Rabbit

相关类别: 一抗

反应种属: Human

标 记 物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	Implicated in the regulation of circadian rhythms through a utocrine and/or paracrine actions By similarity.
Applications:	WB
Name of antibody:	NMS
Immunogen:	Synthesized peptide derived from internal of human NMS.
Full name:	neuromedin S
Synonyms :	0
SwissProt:	Q5H8A3
WB Predicted band size:	18 kDa
WB Positive control:	MCF-7 cells, 293 cells, Huvec cells, HepG2 cells and Jurkat cells lysates
WB Recommended dilution:	500-3000

