

## 兔抗 NMU 多克隆抗体

中文名称：兔抗 NMU 多克隆抗体

英文名称： Anti-NMU rabbit polyclonal antibody

别名： neuromedin U

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NMU

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

<b>Background:</b>	This gene encodes a member of the neuromedin family of neuropeptides. The encoded protein is a precursor that is proteolytically processed to generate a biologically active neuropeptide that plays a role in pain, stress, immune-mediated inflammatory diseases and feeding regulation. Increased expression of this gene was observed in renal, pancreatic and lung cancers. Alternative splicing results in multiple transcript variants encoding different isoforms. Some of these isoforms may undergo similar processing to generate the mature peptide.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	NMU
<b>Immunogen:</b>	Fusion protein of human NMU
<b>Full name:</b>	neuromedin U
<b>SwissProt:</b>	P48645
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	100-300

